

Reynard

Premier Disinfectant Wipes



Premier Disinfectant Wipes

Reynard's Premier Disinfectant Wipe has a white apertured non-woven substrate. This allows for extremely good traction and pick up of dirt/bacteria. Most other surface wipes are made from thermal bond polyester which push the dirt/bacteria around, rather than pick it up.



Directions: Use the required number of wipes.
Dispose of used wipes properly; do not flush.



Product specifications and availability may vary according to country.
Always read the label and product information before use.

Supporting documentation: Full technical specifications are available upon request.
Bespoke solutions: Should your market require more or fewer wipes, lightly scented alternatives or even biodegradable wipes, for example. Let us know.
We are more than happy to develop bespoke solutions.

1 step cleaning & disinfection

Hygienic, easy to use, convenient & efficient. Reynard Premier Disinfectant Wipes will clean and disinfect most hard surfaces and equipment in one easy step. They are safe to use on mattress covers, beds, tables, work tops, trays, trolleys and other hospital equipment. **MPI C44 & Howard Wright Approved**

Fast and effective

Reynard Premier Disinfectant Wipes will kill 99.9% of common germs and have been tested to be effective against a wide range of bacteria. They are

ALCOHOL, FRAGRANCE & CHLORINE FREE .

Cost-effective

Reynard Premier Disinfectant Wipes will clean and disinfect in one. There is no need for extra products. Reynard wipes eliminates the need for rags, spray bottles and chemicals, thus saving on laundering, chemical and storage costs. The non tainting, low foaming formula does not need to be rinsed from surfaces and the large wipe size of 33 x 22cm ensures you use less wipes for all applications.

Technical Performance:

Bactericidal—MRSA, VRE, E Coli, Pseudomonas Aeruginosa, Staphylococcus aureus, Listeria, Salmonella & Enterococcus hirae,

Mycobactericidal—Mycobacterium avium, Mycobacterium terrae,

Fungicidal—Aspergillus Niger,

Virucidal—Influenza A (H1N1), HIV, Hep B, Vaccina, Norovirus,

Sporicidal—C Diff

This Multi Surface Disinfection product has been tested as per EU Standards, EN1276, EN13727, EN14348, EN13704, EN1650, EN14675, EN14476, EN14562 and EN14563

Reynard Premier Disinfectant Wipes



One step for cleaning and disinfecting. Reynard Premier Disinfectant Wipes are fragrance, Alcohol and Chlorine Free to avoid any risk of tainting. Economical – disposable, ready-soaped, rinse-free and large wipe size brings an end to bowls of water, detergents and cloths. Being portable makes cleaning and disinfecting a lot easier.



RHS 216 Premier Disinfectant Wipes

Technical specifications

- 33 x 22cm
- 50 wipes per pack, 12 packs per carton

Dispenser

- Clear perspex dispenser



RHS 218 Premier Disinfectant Wipes

Technical specifications

- 33 x 22cm
- 100 wipes per pack, 6 packs per carton

Dispenser

- Clear perspex dispenser



FORMULATION:

Plantacare 1200 UP Biodegradable non-ionic surfactant with excellent non-streaking and superior wetting properties to provide the wipe with an effective cleaning action and enable the liquid to effectively wet out surfaces.

Microcare IP20 Fungicide to protect packs of the wipes against fungal spoilage due to contamination being introduced whilst a pack of the wipes is being used.

Tetradodium EDTA Chelating agent to enhance the biocidal activity of the wipe.

Bardac The active ingredient (4%). The biocidal components derived from 4% Bardac 114 are as follows:

Didecyl dimethyl ammonium chloride

C12/18 alkyl dimethylbenzyl ammonium chloride

C12/14 alkyl dimethylethylbenzyl ammonium chloride

Irradiated deionised water

The recommendations regarding our preparations are based on scientific tests and are given in good faith. More detailed recommendations, e.g. Regarding material compatibility are only possible in particular cases. Our recommendations are without obligation and do not constitute a warranty. They do not preclude the company's own testing for the intended purposes and processes. In this respect we cannot accept any liability.

This complies with our general conditions of sale and supply.

United Kingdom

Boston House
Grove Technology Park
Downsview Road, Wantage
Oxfordshire OX12 9FF
+44 845 263 7335

Australia

Suite 403, 44 Hampden Road
Artarmon
Sydney, NSW 2064
enquiries @reynardhealth.com
☎ 1300 667 199

New Zealand

17 Napier Road
Havelock North
Hawke's Bay
4130 New Zealand
☎ +64 (0)6 650 0709

International

17 Napier Road
Havelock North
Hawke's Bay
4130 New Zealand
+64 (0)6 650 0709

www.reynardhealth.com